

2X CTAB DNA Isolation Protocol to separate Host and Symbiont Tissue
(Modified from Coffroth et al 1992)

PREPARATION:

1. Make sure 65° C waterbath is turned on and water level adequate on day of extraction. Put 2X CTAB in H₂O bath to thaw.
2. **If samples are frozen:**
 - a. Label eppies with sample numbers, if necessary (label one set as "Z" and one as "A")
 - b. Place eppies in a rack and add 300 ul of 2x CTAB buffer
 - c. Chill two pair of forceps, two razor blades and a metal "cutting board" wrapped in aluminum foil.
 - d. Keeping samples on dry ice, remove coral piece from its cryovial with a pair of well-chilled forceps. **BE VERY CAREFUL TO PREVENT THE SAMPLE FROM THAWING AT THIS TIME.** Switch to the other pair of forceps if necessary. Scrap off the of sample and **place it in appropriate "Z" eppie.** Wipe each instrument used and the cutting surface with a Kimwipe in between samples to prevent cross-contamination. Keep tissue on dry ice at all times; store at -70°C
3. **If extracting samples stored in ETOH,** transfer a portion of the tissue to 1.5ml eppie with 300 ul 2X CTAB buffer and proceed with Protocol below.

PROCEDURE:

- 1) Add 300µl 2X CTAB to the eppie labeled "Z" prior to adding the tissue and grind sample with blue dounce, rinse dounce w/ 300µl 2X CTAB into eppie (total CTAB volume 600µl.)
- 2) Place the "Z" eppie in a 15 ml round bottom Sarsdarst tube and spin in the floor (low speed) centrifuge at 500 rpm for 5 min.
- 3) Remove the supernatant from the "Z" and place in the appropriate "A" tube;
- 4) Add 300 ul of 2X CTAB to the "Z" tube, vortex briefly and repeat step (2)
- 5) Discard the supernatant and repeat step (4)
- 6) Discard the supernatant and add 300 ul of 2X CTAB to the "Z" tube.
- 7) Repeat step (5) and (6)
- 8) After above spin, resuspend the pellet in 300 ul 2X CTAB
- 9) Add one layer of glass beads to the cells and rupture cells by placing on Turbo-mixer (vortex mixer) for 5 minutes. .
- 10) Add 3.0µl of Proteinase K (stock=20mg/ml) to the tubes labeled "A" and 1.5 µl of Proteinase K to the tubes labeled "Z". Mix well by pipeting up and down, **don't vortex.** Incubate at 65°C for 30 min. Record time in and out of water bath. Invert tubes intermittently if possible.
- 11) Add equal or greater volume of chloroform-isoamylalcohol (CIA - 24:1), (~600ul for the "A" tubes and ~300ul for the "Z"). Buzz each sample on vortex machine

- for approximately 15 seconds. Centrifuge for 5 minutes at 12K (high speed in microfuge).
- 12) Transfer the aqueous phase (top) with P200 to labeled 1.5 ml clear eppies.
 - 13) Add an equal or greater volume of buffered phenol: chloroform: isoamylalcohol (ϕ : CIA - 25:24:1) (~600ul for the "A" tubes and ~300ul for the "Z"). Buzz each sample on vortex machine for approximately 15 seconds and centrifuge for 5 minutes at 12K.
 - 14) Transfer aqueous phase (top) to final tubes. Repeat the phenol/chloroform step if the interphase is large or the aqueous layer is not clear.
 - 15) Add 2 volumes (1000ul for the "A" tubes and ~500ul for the "Z") of cold 95% ETOH, shake it very well by inverting and place at -20°C to precipitate the DNA over night or at -70°C in hall for ½ hour.
 - 16) Collect precipitate by centrifugation (high speed for 30 minutes). NOTE: In the case of zoox from polyps, often you cannot see a pellet. It is essential that you always orient the tubes in the same way so that you know where the DNA is. Discard supernatant, but be careful not to lose pellet – pipet or pour out supernatant from opposite side from where DNA should be.
 - 17) Add 500ul 70% ETOH, mix and spin at 10k for 5 minutes. Remove the ETOH being carefully not to lose the pellet.
 - 18) Repeat the 70% ETOH wash and spin for 5 minutes. Again remove ETOH, spin the samples for 30 sec and remove the rest of the ETOH with a pipet. Use a greater volume of ETOH if pellet is large or initial ETOH ppt. volume greater.
 - 19) Dry the pellet in freeze drier or speed-vac for 15 minutes.
 - 20) Resuspend in 5 to 15ul TE, depending on size of pellet. To do this, carefully wash the side of the tube where the DNA should be repeatedly with the 5 to 15ul of TE that you are resuspending in. The idea is to dissolve all DNA adhering to the side of the tube in a minimum volume. (If there is no pellet, then 10ul. If there is pellet, then 15ul). If it doesn't dissolve, incubate at 65°C for 5 minutes, vortex briefly and spin down.
 - 21) Store samples at 4°C freezer (short-term, i.e., overnight in well-labeled red bullet boxes) or -20° C long-term in white boxes.